

Digital Inspection video Microscope JBM-MZ50





Digital Monozoom microscopes adopts apochromatic parallel optical imaging system and provide high resolution and sharp stereoscopic images, 4K HDMI digital camera with measurement function can work without PC.

This microscope is designed for the applications of circuit board, electronics, semi-conductor and related industrial inspection and scientific research areas. Modularization product portfolio and excellent optical performance make it a better choice in these areas.

Technical Specifications

Brand Name	JB MICROSCOPE
MODEL	JBM-MZ50
Zoom Ratio	01 :07
Zoom objective	0.7X-5X
Eyepiece	0.4x
Working Distance	100mm
Base Vertical post height	240mm
Vertical post diameter	φ32mm
Microscope body holder	φ50mm
Focus travel	50mm
Base size	320x260x16mm
Material	Metal



+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++