

Digital Microscope JBM-AD10i





HDMI Output
UHD 1920×1080 24FPS

USB Output
PC software with
Measurement



Andonstar's Innovative Pro Boom Arm Stand



Specification:

Model Number	AD10i
Imaging sensor:	4 mega pixels high definition sensor
screen size	10 inch
Magnification ratio:	up to 2000X (high-definition multimedia interface monitor 27 inch)
Video output:	UHD 2880*2160 24FPS; FH-D 1920*1080 60 FPS/30FPS, H-D 1280*720 120FPS
Photo resolution:	Max. 4032*3024 (12M 4032*3024)
Frame rate:	Max. 120f/s (under 600 lux brightness&H-DP120)
PC support:	YES
Photo format:	JPG
Minimum focus distance:	Micro-SD card, up to 32G (not included)
Minimum focus distance:	4mm
Power source:	5V DC
Upgraded Stand size:	20cm*18cm*33cm

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++