

Digital portable Inspection Microscope JBM-DM10



HDMI Digital Microscope with 10.1" IPS Screen

Capture every moment
Share every experience



FHD 1080P IPS Screen



10.1"
IPS Screen

HD

FHD
Full HD Display



178°
Viewing Angles

Make the Best of
Everything



Video
1080P



Photo
20MP max



Audio
Lossless audio sound

Remote Control


Get Rid of the "Line" Control



1m-5m
Remote control




No more shake



Technical Details

Model Number	JBM-DM10
Theory	Video Microscope
Drawtube	Monocular
Material	ABS
Dimension	31*29*10 Cm
Magnification	2000x
Type	Digital Microscope, Stereo Microscope, Electron Microscope
Screen	10.1 inches IPS Display
Display Resolution	1024*600 Resolution
Video Maximum Pixel	1080p
Photo Resolution	20M,16M,12M,8M,5M,2M
Video Format	MP4
PC Resolution and Image Transfer Rate	1280*720/25fps
PC Operating System	Win7, Win 8.1, Win10,Mac Os X 12 or higher
Interface and Signal Transmission mode	Micro/USB3.0/HDMI
Polymer Lithium Battery	103665-3000mah



Multiple Output Modes

① Use alone




② Memory card



③ USB output



④ HDMI output



+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++