

Digital portable Inspection Microscope JBM-DM12



High-Definition Multimedia Interface



10.1" IPS Screen

178° full viewing angle clear



10.1" IPS Screen

178° View Angle

90° Screen Folding

Compatible with Multiple Systems



Real-time Voice Recording While Exploring


Increased productivity

Ease of sharing



Technical Details

Model Number	JBM-DM12
Theory	Video Microscope
Drawtube	Monocular
Material	ABS
Dimension	31*29*10 Cm
Magnification	2000x
Type	Digital Microscope, Stereo Microscope, Electron Microscope
Screen	10.1 inches IPS Display
Display Resolution	1024*600 Resolution
Video Maximum Pixel	1080p
Photo Resolution	20M,16M,12M,8M,5M,2M
Video Format	MP4
PC Resolution and Image Transfer Rate	1280*720/25fps
PC Operating System	Win7, Win 8.1, Win10,Mac Os X 12 or higher
Interface and Signal Transmission mode	Micro/USB3.0/HDMI
Polymer Lithium Battery	103665-3000mah



++++
++++
++++
++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

++++
++++
++++
++++