


Stereo Zoom Microscope

JBM-750B



Stereo zoom microscopes allow quick zoom from low to high magnification. Using a stereo zoom microscope makes it simple to view a sample in the entire field of view, then focus and zoom in on a particular part of interest .



[illegible]

Brand Name	JB MICROSCOPE
MODEL	JBM-750B
Eyepiece	WF10X/22mm
Zoom Objective	0.7-50X
Zoom Range	7.1 : 1
Working Distance	100mm
Head Mounting Size	76mm in diameter
Dust cover	Yes (Included)
Vertical post height	240mm
Vertical post diameter	φ32mm
Microscope body holder	φ50mm
Focus travel	50mm
Vertical post height	240mm
Base size	320x260x16mm
Material	Metal

++++
++++
++++
++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

++++
++++
++++
++++