


Auto focus Digital portable Inspection Microscope JBM-A05+



JBM-A05+ Auto Focus 1- 500X ZOOM HD 5MP CMOS SENSOR USB
MINI ELECTRON POCKET VIDEO MICROSCOPE CAMERA PRICE WITH STAND



- Theory: Video Microscope
- Drawtube: Monocular
- Type: Pocket Microscope, Digital Microscope
- Customized Support: OEM, ODM, OBM
- Model Number: A005+
- Material: Aluminum
- Dimension:18x12x4.4 cm
- Magnification:1-500X
- Focus: Auto Fine Focus & Manual Focus
- Magnification: 1- 500x (Continuous Zoom)
- Sensor: HD 5 MP
- Resolution:2592×1944
- Frame Speed: Max30 Frames/ Second
- Port: USB 2.0
- Take photo: yes
- Record Video: yes
- Calibration: yes



+++++

+++++

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

WWW.JBM ICROSCOPE.COM

INFO@JBINSTRUMENTS.COM

+91 9354741245

+++++