



Microscope Adapter 0.5X C-Mount Reduction Lens

Microscope Adapter, 0.5X C-Mount Reduction Lens Eyepiece Lens 23.2mm Mounting For Camera Digital Eyepiece

Features:

0.5X magnification, which makes it easier to observe small objects more clearly.

C-mount adapter lens, mainly used to camera or digital eyepiece to the microscope.

This product is glass in use and long service life.

Well-made, and use, can provide clear images- for easy observation.

microscope accessories, excellent make sure better observation effect.



Large image but small view




Image is smaller, view is larger



Specifications:

Metal

Mounting Size: 23.2mm (to binocular mount)

C-mount Thread: 25.4mm (to electronic eyepiece or camera)

Suit Microscope: biological microscope and stereo-microscope

Suitable 1/2" 1/2.5"

