


DIGITAL MICROSCOPE CAMERA

JBM-1000 (USB3.0)

8MP | High Resolution
Crystal Clear Imaging



8MP
High Resolution



USB 3.0
High Speed



Ultra HD
Image Quality



Wide
Compatibility



Sharp & Clear
Detail



Advanced CMOS
Sensor



Plug & Play
Setup



Multi-Platform
Support



USB3.0 high speed transmission interface, simple and convenient at the same time to ensure stable high transmission rate full resolution can realize 30 frames/s fast imaging, up to 60 frames/s



8 Mega
Pixel



USB 3.0
High Speed



High Performance
Sensor



Crystal Clear
Imaging



Up to 60fps
High Frame Rate

TECHNICAL SPECIFICATION

Brand	JB MICROSCOPE
Model number	JBM-1000 (USB3.0) 8MP
Camera Pixel	8 Mega Pixel
Chip Model	IMX678
Optical Size	1/1.8 inch
Image size	2.0μm X 2.0μm
Resolution	3840(H) X 2160(V)
Camera frame rate	up to 60 fps
Output Interface	USB3.0 high speed output
External Device	Desktop PC or PC Notebook
Optical Interface	Standard C connector
Auto Settings	Exposure, Brightness, White Balance
Manual Settings	Exposure, Gain, Noise Reduction, Gamma



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245