

LCD Microscope Digital Camera JBM-L7





- New generation 7-inch large display with a built-in camera, built-in wi-fi transmitter.
- The mobile phone can be connected to wi-fi, synchronous observation.
- 7-inch high-definition large screen, the image is clearer.
- You can rotate it by 90 degrees
- The ultra-long working distance is conducive to spot welding by mobile phone maintenance, circuit board inspection, etc.
- Remote button remote control/External fill light/Free reflector.
- Continuous zoom trinocular simultaneous observation, clear image/easy operation.
- Large wide-angle lens/continuous zoom/adjustable LED/diopter adjustment

Brand	JB Microscope
Screen Size	7 Inch HD LCD Display
Resolution	1080P
Brightness	Adjustable
Interface	Micro/USB 2.0, HDMI
Power Supply	3000mAH
Memory Card	16 GB TF Card Include (Max support 128G card)
Output Port	HDMI
Viewing angle	IPS full view
WIFI distance	Ten meters



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com | info@jbinstruments.com | +91 9354741245