



LCD Microscope Camera JBM-L9





Technical Parameters

Brand	JB Microscope
Model no	JBM-L9
Series	Integrated Display camera
Sensor Size	1/2.5" CMOS
Pixel Size	2.2*2.2um
Effective pixels	5 million pixels
Resolution	2048*1536@15FPS(5MP)
Screen Resolution	2048*1536 Retinal Screen
Operating system	Built-in Android System
Output	USB2.0 TF card
Output Ports	1x USB 2.0 Port
Power	DC-12V/2A



[illegible]

- **UI:** Easy mouse control (Optional Bluetooth mouse), Android5.1 operation system
- **Mark function:** Point coordinates, Reticle, Coordinates, Text annotation
- **Length measuring:** Straight length, broken line length, curve line length, the distance between parallel lines, point-line distance
- **Center distance of circles:** Center distance of circles
- **Geometry measurement:** Line length, radius setting circle, two points setting circle, three points setting circle, concentric circles e, concentric circles
- **Geometry area:** Polygon and square
- **Image capture:** Support 16 million static image capture
- **Video record:** Support 1080P video
- **Measuring tool:** S-EYE.

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++