


LCD Microscope Digital Camera

JBM-L10





Technical Details

Brand	JB MICROSCOPE
Model number	JBM-L10
Screen Size:	10 Inch HD LCD display
Resolution	1080P
Brightness	Adjustable
Interface	Micro/USB 2.0, HDMI
Power supply	3000mAh
Memory card	16 GB TF Card Include(Max support 128G card)
Output port	HDMI
Viewing angle	IPS full view
Wifi Distance	Ten Meter



+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++