




JB INSTRUMENTS PVT. LTD.
"Delivering Innovative Solutions To Life's Everyday Needs"

Textile Soft

Textile Soft



Specifications

- This Software useful for micro measurement of fiber, yarn, threads, textile material, non-woven fabrics etc.
- The real full screen display of image, acquisition, annotation, enhancement with hundred of filters, measurements, statistical processing of measured data are few of the capability of Textile Soft.
- Photographs of all fibers view are pre-loaded in the software to identify fibers without prior knowledge in the same. The library can be created & retrieved within the same program.
- Micro analysis of yarns and fiber sections.
- Measurements of length, surfaces, perimeters, angels, distance between two points or lines or objects.

JB INSTRUMENTS PVT. LTD.

For more information contact us www.jbmicroscope.com

J-9/36, Unit No. 302, First Floor,
Rajouri Garden, New Delhi-110027

Cell : +91-9873571533, +91-9354741245
E-mail : info@jbinstruments.com

www.jbinstruments.com, www.jbmicroscope.com



JB INSTRUMENTS PVT. LTD.
"Delivering Innovative Solutions To Life's Everyday Needs"

- **Fast and easy way the fitness analysis of single fiber.**
- **Check the purchased material can be identify the type of fiber, comparing it with the fiber pictures stored in the data bank with the longitudinal and the sectional views.**
- **Check and measures the quality and shape of Lycra or synthetic multifilament single threads.**
- **Analyze the compactness of non-woven fabrics.**
- **Analyze the Yarn structure and detect possible defects.**
- **Detect, identify and measure possible impurities contained in textile materials.**
- **Analyze the compactness of no-woven fabrics.**
- **Measure section surface and perimeters.**
- **Analyze mechanical part like needle points spinnerets etc.**
- **Macro analysis of a yarn section.**
- **Reduce the fabric warp and weft density to a cm or inch.**

JB INSTRUMENTS PVT. LTD.

For more information contact us www.jbmicroscope.com

**J-9/36, Unit No. 302, First Floor,
Rajouri Garden, New Delhi-110027**

**Cell : +91-9873571533, +91-9354741245
E-mail : info@jbinstruments.com**